



Serve

from  **psl**

Spares, Service and
Customer Care

Operators Require Downtime. **Your Production Does Not.**



Your bellows and agitator assembly are critical for your production.

All bellows have a dedicated cycle life and as such, must be replaced to prevent loss of a batch.

If your bellows has exceeded its cycle life, it can take up to 16 weeks to replace and your equipment cannot run without them. These bespoke consumable items are manufactured to order.

PSL advises our clients to hold items for an agitator and bellows assembly in stock at your facility so this maintenance can be planned in advance.

**Not having on-site
spares WILL cost you
weeks of lost production.
Can you afford this?**



Please contact Serve@powdersystems.com or call
+44 (0) 151 448 7700 for more information and a quotation.



Serve

from  **psl**

Spares, Service and
Customer Care

Preventive Maintenance for Reliable Manufacturing

A Preventive Maintenance programme ensures equipment performance and minimises production downtime by taking care of the equipment before problems arise.

Preventative maintenance can be thought of as being like taking your car for a service; we look after each piece of equipment, conducting thorough inspections and if needed replacing components.

Take Control

Our dedicated service engineers will accommodate your production schedules, including shut down periods, enabling you to take control of maintenance activities when convenient to you.

Planned Maintenance

PSL's Preventive Maintenance programme is a simple

method to ensure that you get reliable and efficient performance from your plant and equipment. It is important to maintain machinery at optimal performance so you don't suffer a critical break down when you need it most.

Our Maintenance Packages range from our Standard Level, which comprises the provision of labour and equipment inspections, through to our Comprehensive Level including all consumable parts and labour for a trouble-free maintenance programme.

- **Optimal Performance**
- **Minimize Interruption**
- **Increase Productivity**

The PSL Serve team can upgrade and refurbish any filtration, drying or high containment installation around the world. Contact us to find out more.



Our dedicated
customer care experts

Technical Support

Spares

Re-qualification
and Upgrade

IQ/QQ Equipment
Packages

Containment
Performance
(OEL) Testing

Glove Integrity
Testing System

On-Site Services

Planned Maintenance

Emergency
Breakdown Cover

Training

Please contact Serve@powdersystems.com or call
+44 (0) 151 448 7700 for more information and a quotation.